



# EVOTEC

## Furniture Catalogue

 **ceasa**  
THE CLINICAL ENGINEERING ASSOCIATION  
OF SOUTH AFRICA  
corporate member

658 Alouette Street  
Elarduspark  
Pretoria

P.O. Box 4333  
Rietvalleirand  
0181

Tel : 012 345 1176 / 4738 / 6421  
Fax : 012 345 4381 or 086 550 4126  
E-mail: [sales@evotec.co.za](mailto:sales@evotec.co.za)

[www.evotec.co.za](http://www.evotec.co.za)

## THE FACTORY



Pictures are for illustrative purposes only and may vary from actual product.

**Dressing Trolley  
Medium - Stainless****Part No: EVO 53/SS/610**

- Stainless Steel trolley
- Stainless Steel shelves
- Guard rail on 3 sides
- 4 x 75mm castors

**Dressing Trolley  
Small - Stainless****Part No: EVO 53/SS/457**

- Stainless Steel trolley
- Stainless Steel Shelves
- Guard rail on 3 sides
- 4 x 75mm castors

**Dressing Trolley  
Large - Stainless****Part No: EVO 53/SS/915**

- Stainless Steel trolley
- Stainless Steel Shelves
- Guard rail on 3 sides
- 4 x 75mm castors

**Dressing Trolley  
Medium - Epoxy****Part No: EVO 53/EP/610**

- Epoxy Coated trolley
- Stainless Steel shelves
- 4 x 75mm castors

**Dressing Trolley  
Small - Epoxy****Part No: EVO 53/EP/457**

- Epoxy coated trolley
- Stainless Steel Shelves
- 4 x 75mm castors

**Dressing Trolley  
Large - Epoxy****Part No: EVO 53/EP/915**

- Stainless Steel trolley
- Stainless Steel Shelves
- 4 x 75mm castors

**Instrument Trolley  
Stainless****Part No: EVO 53/EP/457**

- Stainless Steel trolley
- Stainless Steel Shelves
- 4 x 75mm swivel castors

**Dressing Trolley  
2 Drawer****Part No: EVO 54**

- Stainless Steel trolley
- Stainless Steel Shelves
- 4 x 75mm swivel castors

**Anaesthetic Trolley  
1- Drawer****Part No: EVO 55**

- Stainless Steel shelf
- Pull out drawer
- Chrome plated Guard rail

**Linnen or Food trolley  
3 Tier****Part No: EVO 63 EP**

- 940 x 900 x 530 mm
- Mild Steel- Epoxy
- 3 Tier mild steel epoxy coated trays
- 4 x 75mm swivel castors

Pictures are for illustrative purposes only and may vary from actual product.



## Two Tier Trolley

Part No: EVO 61 EP

- Mild steel Epoxy coated
- 2 Stainless Steel trays
- 4 x 75mm swivel castors



## Tea Trolley

Part No: EVO 115

- Specialized Item - Made on request



## ICU Bed End Trolley

Part No: EVO 105

- Mild steel Epoxy coated
- Formica wood top with aluminium edging
- Top drawer ( optional 4 drawer left with and door on right
- Needle box holder , soap dispenser , x-ray holder



## Emergency Trolley 5 Drawer

Part No: EVO 26E

- Mild steel Epoxy coated
- Stainless Steel top with rail
- 5 lockable drawers
- Drip stand holder , oxygen holder, CPR board
- 4 x 125mm castors (2x brake)



## Anaesthetic Trolley 6 Drawer & 3 Shelves

Part No: EVO 26C

- Mild steel Epoxy coated
- Stainless Steel top
- 6 drawers & 3 shelves
- Drip stand holder , oxygen holder, CPR board
- 4 x 125mm castors



## Anaesthetic Trolley 6 Drawers

Part No: EVO 26C/A

- Mild steel Epoxy coated
- Stainless Steel top
- 6 drawers
- 4 x 125mm castors



## Drug Cabinet

Part No: EVO 26 B/610

- Mild steel Epoxy coated
- 2 x Lockable doors
- Writing shelf
- Stainless steel tray with bowl holders. Stainless tray with rail
- Castors



## Neo-Natal Bed End Trolley

Part No: EVO 106

- Mild steel Epoxy coated
- Formica wood
- 1 x Drawer
- 1 x Wire basket
- 4 x 75mm castors



## Child Cot

Part No: EVO 7

- Mild steel Epoxy coated
- Drop down sides
- 4 x Castors with 2 x brakes
- Mattress sold separate



## Bassinette

Part No: EVO 19

- 700L x 440W x 820H mm
- Clear acrylic crib
- Water resistant taffeta cover mattress
- 4 x swivel castors

Pictures are for illustrative purposes only and may vary from actual product.



### Overbed Table Fixed Height Stainless

Part No: EVO 22/C

- Stainless Steel frame & top
- L 1090 x W 400 x H 900mm
- 4 x 75 mm swivel castors



### Overbed Table Fixed Height Wooden top

Part No: EVO 22/A

- Mild steel epoxy coated
- L 1090 x W 400 x H 900mm
- 4 x 75 mm swivel castors
- Solid formica top



### Overbed Table Fixed Height Epoxy

Part No: EVO 22

- Mild steel epoxy frame & top
- L 1090 x W 400 x H 900mm
- 4 x 75 mm swivel castors



### Hi-Lo Overbed Table Epoxy

Part No: EVO 22 D

- Mild steel epoxy frame & formica top
- Gas spring adjustment to height between 800 to 1100 mm
- 4 x 75 mm swivel castors



### Hi-Lo Overbed Table Stainless Steel Top

Part No: EVO 22 DS

- Mild steel epoxy frame & Stainless Steel Top
- Gas spring adjustment to height between 800 to 1100 mm



### Locker

Part No: EVO 30 A

- Mild steel epoxy & Formica wood top
- Locker with Drawer
- 4 x swivel castors



### Locker- Single

Part No: EVO 28 A

- Mild steel epoxy & Formica wood top
- Locker with 1x shelf
- 1 x Drawer
- Epoxy coated towel rail
- 4 x swivel castors



### Locker- Single

Part No: EVO 28

- Mild steel epoxy & Stainless Steel top
- Locker with 1x shelf
- 1 x Drawer
- Epoxy coated towel rail



### Combo Locker- Short

Part No: EVO 109

- Mild steel epoxy & Formica top
- Cupboard with roller catch and shelf
- 1 x Drawer
- Epoxy coated towel rail



### Combo Locker- Long

Part No: EVO 109L

- Mild steel epoxy & Formica top
- Cupboard with roller catch and shelf
- 1 x Drawer
- Epoxy coated towel rail

Pictures are for illustrative purposes only and may vary from actual product.



## Mayo Table with Large Tray

Part No: EVO 62

- Stainless Steel frame with stainless tray
- Height adjusted 800 to 1300 mm
- 4 x 50mm swivel castors



## Bed Screen

Part No: EVO 74

- Mild steel epoxy coated
- Panel size 1700 x 600mm
- Available in 3, 4, 5 panels
- Curtain sold separate



## Double Bowl Stand

Part No: EVO 33

- Stainless Steel
- H835 x L610 x W390 mm
- Bowl Size rings 295mm
- 4 x swivel castors
- Bowls sold separate



## Single Bowl Stand

Part No: EVO 34

- Stainless Steel
- H835 x L390 x W390 mm
- Bowl Size ring 295mm
- 4 x swivel castors
- Bowl sold separate



## Kick-About Bucket type

Part No: EVO 35 EP

- Epoxy coated mild steel
- Single bucket ring
- 3 x swivel castors
- Bucket sold separate



## Kick-About Bucket type

Part No: EVO 35 SS

- Stainless Steel
- Single bucket ring
- 3 x swivel castors
- Bucket sold separate



## Kick-About Bowl type

Part No: EVO 36 SS

- Stainless Steel
- Single bowl ring
- 3 x swivel castors
- Bowl sold separate



## Drip Stand-Mobile Stainless Steel

Part No: EVO 68 SS

- Stainless Steel
- 5 x swivel castors
- Adjustable height 1300 to 2000 mm



## Drip Stand-Mobile Epoxy coated

Part No: EVO 68 EP

- Epoxy Coated
- 5 x swivel castors
- Adjustable height 1300 to 2000 mm



## Drip Stand-Mobile Heavy Base

Part No: EVO 68 HB

- Manufactured to customers specifications
- 5 x swivel castors

Pictures are for illustrative purposes only and may vary from actual product.



## Oxygen Cylinder Trolley

Part No: EVO 69A

- Epoxy coated mild steel
- 10mm base plate
- Front caster 150mm back castors 75mm
- For 180mm diameter cylinders



## Oxygen Cylinder Stand

Part No: EVO 70

- Epoxy coated mild steel
- Size H680 x L400 x W400mm
- For 180mm diameter cylinders



## Linen Receiver

Part No: EVO 76

- Epoxy coated mild steel
- Size W610 x H900mm
- 610 mm diameter ring
- Swivel castors
- Canvas bag sold separate



## Linen Receiver with 3 shelves

Part No: EVO 77

- Epoxy coated mild steel and shelves
- 445mm diameter ring
- Swivel castors
- Canvas bag sold separate



## Clinibin 85L

Part No: EVO 99

- Epoxy coated mild steel
- Size H1080 x W530 x L520mm
- Swivel castors
- Removable plastic bin with handle



## Clinibin 26L

Part No: EVO 99 A

- Epoxy coated mild steel
- Size H890 x W490 x L520mm
- Swivel castors
- Removable plastic bin
- Bin lid activated by pedal



## Monitor trolley

Part No: EVO Cust

- Made on request to customers specifications



## Anaesthetic Stool

Part No: EVO 78

- Epoxy coated mild steel
- Aluminium seat
- Swivel castors 4 x 50mm



## Operating stool

Part No: EVO 78 A

- Epoxy coated mild steel
- Adjustable upholstered seat top
- fixed



## Operating stool

Part No: EVO 78 B

- Epoxy coated mild steel
- Adjustable upholstered seat top
- Swivel castors 4 x 50mm

Pictures are for illustrative purposes only and may vary from actual product.

**Visitors Bench****Part No: EVO 79**

- Epoxy coated mild steel
- 3 sizes 915 , 610, 450 mm
- 20mm laminated wood seat
- 30mm x 100mm carry slot

**Foot Stool Single****Part No: EVO 80**

- Epoxy coated mild steel
- Non slip rubber top
- Aluminium edging

**Foot Stool Double****Part No: EVO 80 A**

- Epoxy coated mild steel
- Non slip rubber top
- Aluminium edging

**Foot Stool Single****Part No: EVO 80 SS**

- Stainless Steel
- Aluminium with pattern tread plate top
- Aluminium edging

**Foot Stool Double****Part No: EVO 80 A SS**

- Stainless Steel
- Aluminium with pattern tread plate top
- Aluminium edging

**Patient Trolley with Stainless Steel Top****Part No: EVO 86**

- Epoxy coated frame
- Stainless Steel top and handles
- Drip stand, oxygen cylinder holder, document holder, 2 brake and 2 swivel castors

**Bed - 2 section Hi-Lo bed****Part No: EVO 1700**

- Two section h-lo bed
- Hydraulic foot pedal mechanism
- Removable Head and foot boards
- 2 x Brake and 2 X swivel castors

**Examination Couch****Part No: EVO 88**

- Epoxy coated mild steel
- Ratchet backrest
- Vinyl covered top

**Examination Couch****Part No: EVO 88 A**

- Epoxy coated mild steel
- Ratchet backrest
- Vinyl covered top
- Breathing hole

**Foldable Couch****Part No: EVO 110**

- Made on request

Pictures are for illustrative purposes only and may vary from actual product.



### Hostel Bed

Part No: EVO 111

- Made on request



### Instrument Cabinet

Part No: EVO 40

- Epoxy coated mild steel
- Lockable
- Double door
- 2 x Adjustable glass shelves



### Poison Cabinet Large

Part No: EVO 45 A

- Epoxy coated
- Lockable door
- 1 x shelf



### Poison Cabinet

Part No: EVO 45

- Epoxy coated
- Lockable door
- 1 x shelf



### Drug cabinet Double door

Part No: EVO 45 E

- Epoxy powder coated
- 1 x Fixed shelf
- 1 x adjustable shelf
- Lockable
- Schedule 7 lockable cabinet inside



### Poison Cabinet Double (on Top)

Part No: EVO 45 D

- Epoxy coated
- 2 x Lockable door
- 



### Deadly Drug cabinet Small

Part No: EVO 45F

- Epoxy powder coated
- Lockable door inside



### Deadly Drug cabinet Large

Part No: EVO 45AA

- Epoxy powder coated
- Lockable door inside
- Shelf inside
- Schedule 7 lockable cabinet inside



### X-Ray Viewing Box

Part No: EVO 90

- Epoxy powder coated
- X-ray viewer



### Chart Holder

Part No: EVO 92

- Epoxy Powder coated
- 2 pockets
- Double hook with protective plastic


Pictures are for illustrative purposes only and may vary from actual product.

**Chart Holder  
Stainless Steel****Part No: EVO 92 SS**


- Stainless Steel
- 2 pockets
- Double hook with protective plastic

**Lifting Pole  
(Monkey pole with handle)****Part No: EVO 20**

- Epoxy coated mild steel frame
- Triangular shaped handle with chain link
- Adjustable

**Urine and Catheter  
bag holder****Part No: EVO 98 EP**

- Epoxy coated mild steel

**Urine and Catheter bag  
holder****Part No: EVO 98 SS**


- Stainless steel

**Cotton Wool  
dispenser****Part No: EVO 48A**

- Epoxy powder coated

**Brush Dispenser****Part No: EVO 48**

- Epoxy coated

**Bed Block****Part No: EVO 15**

- Epoxy powder coated
- To increase bed height

**School Desk****Part No: EVO 111**

- On request only

**Wall Brackets****Part No: EVO 120**

- Epoxy powder coated or Stainless steel
- Made on request to customers specifications

**Please forward requests on all  
items NOT listed**

Pictures are for illustrative purposes only and may vary from actual product.

# EVOTEC

Your one-shop solution  
Since 1997

**DIRECT IMPORTER  
MANUFACTURE AND  
DISTRIBUTOR**

Pictures are for illustrative purposes only and may vary from actual product.

658 Alouette street  
Elarduspark  
Pretoria

P.O. Box 4333  
Rietvalleirand  
0174

Phone: 012 345 1176 / 4738 / 6421  
Fax : 012 345 4381 or 086 550 4126  
E-mail: [sales@evotec.co.za](mailto:sales@evotec.co.za)

[www.evotec.co.za](http://www.evotec.co.za)